




**CARFIL QUALITY**  
YOUR SUPPLIER FOR ALL LABORATORY NEEDS

## Dacron SkinButton: a catheter exit Site Option



### Features

The Dacron® Mesh SkinButton can be used as an alternative to harness or jacket for rodents with externalized catheters. These disposable buttons are perfect for longer term studies as the subcutaneous tissue will grow into the mesh ensuring stability. The catheter passes through the stalk and can be attached to a tether.



**PAVAN SERVICE BVBA**

Beyntellus 3 - B-2360 Oud-Turnhout - Tel. +32 (0)14 45 13 10 - admin@carfil.be - BTW BE 0420.681.872

[www.carfil.be](http://www.carfil.be)