



**CARFIL QUALITY**  
YOUR SUPPLIER FOR ALL LABORATORY NEEDS

## 3-way valves



### Caractéristiques

To be used for simultaneously changing the direction of both the **SUPPLY** of the anaesthetic gas mixture **AND** the **ECHAUST** of the anaesthetic gas mixture. The risk of making mistakes when the gas is supply and gas exhaust have to be changed individually from f.i. the Induction box to the face mask, is reduced substantially bij using these valves.

### 3-Way Valve

By turning the knob from the right to the left, both supply and exhaust are being directed from the face mask to/from the Induction box.

### Numéro de commande

180000159

Double 3-way valve



**PAVAN SERVICE BVBA**

Beyntellus 3 - B-2360 Oud-Turnhout - Tel. +32 (0)14 45 13 10 - admin@carfil.be - BTW BE 0420.681.872

[www.carfil.be](http://www.carfil.be)